

Panasil putty makes further advances.



Panasil® putty

A-silicone precision preliminary
impression material

Panasil putty – our quality putty now has improved properties. With any Panasil putty – regular, fast or soft – frustrating waiting times are a thing of the past. Using the one-step sandwich technique for example, Panasil putty soft and the light bodies, Panasil contact plus or Panasil contact two in one, set simultaneously. Panasil putty retains its tried and tested properties of extreme flexibility, smooth kneading and high elasticity, ensuring that it is still easy to use.



Panasil® putty – Advantages at a glance:

- Excellent setting properties produce short setting times and guarantee a quick impression with an adequate working time.
- Different degrees of flexibility allow it to be used both for the two-step putty wash impression technique (Panasil putty and Panasil putty fast set) and for the sandwich technique (Panasil putty soft).
- The excellent bond between the putty and light body material produces more accurate impressions.
- Its smoothness and easy kneading guarantee effortless handling.

Working and setting times:

Panasil putty:

Total working time: 2 minutes
Intraoral setting time: 2 minutes
Setting time: 4 minutes

Panasil putty fast set:

Total working time: 1 minute 30 seconds
Intraoral setting time: 2 minutes
Setting time: 3 minutes 30 seconds

Panasil putty soft:

Total working time: 2 minutes
Intraoral setting time: 2 minutes
Setting time: 4 minutes

Panasil® putty – Kneadable, addition-curing polyvinyl siloxane precision preliminary impression material.

Panasil putty normal pack, 900 ml, Cat. No. 11101
Panasil putty economy pack, 4 x 900 ml, Cat. No. 11103



Panasil® putty fast set – Kneadable, addition-curing polyvinyl siloxane precision preliminary impression material. Quick-setting.

Panasil putty fast set normal pack, 900 ml, Cat. No. 11141
Panasil putty fast set economy pack, 4 x 900 ml, Cat. No. 11143



Panasil® putty soft – Kneadable, addition-curing polyvinyl siloxane precision preliminary impression material. Low final hardness.

Panasil putty soft normal pack, 900 ml, Cat. No. 11121
Panasil putty soft economy pack, 4 x 900 ml, Cat. No. 11123



A perfect combination

Panasil putty combines perfectly with the light bodies, Panasil contact plus and Panasil contact two in one, to ensure high precision and an extended working time.